**BEYOND BIOTECHNOLOGY  
BEYOND INNOVATION**

FOR BETTER HEALTH CARE

The logo for addbio MEDITEK features the word "addbio" in a bold, lowercase sans-serif font, with a stylized leaf icon above the "i". Below "addbio" is the word "MEDITEK" in a smaller, uppercase sans-serif font.

**addbio**  
MEDITEK

**Research Products Catalog  
for Molecular Biology**

Standard PCR

---

Reverse Transcription

---

Real-time PCR

---

Sample Preparation

2023-2024 Supplement

**Standard PCR Polymerase & Master Mix.** 02

• **Normal Taq DNA Polymerase & Master Mix.**

Add Taq DNA Polymerase  
Add Taq DNA Master (2x conc.)

• **Hot-start Taq DNA Polymerase & Master Mix.**

AddStart Taq DNA Polymerase  
AddStart Taq DNA Master (2x conc.)

• **Multiplex PCR Master Mix.**

AddMulti Taq DNA Master (2x conc.)

• **High Performance PCR Master Mix.**

AddStart HP Taq DNA Polymerase  
AddStart HP Taq DNA Master (2x conc.)

**Reverse Transcription** 07

• **cDNA Synthesis**

AddScript cDNA Synthesis Kit

• **One-step RT-PCR**

AddScript RT-PCR Master (2x conc.)

**Real-time PCR Master Mix.** 10

• **Real-time PCR (Intercalating Dye)**

AddGreen qPCR Master (2x conc.)  
AddGreen RT-qPCR Master (2x conc.)

• **Normal Real-time PCR (TaqMan Probe)**

AddProbe qPCR Master (2x conc.)  
AddProbe RT-qPCR Master (2x conc.)

• **Fast Real-time PCR (TaqMan Probe)**

AddFast qPCR Master (2x conc.)  
AddFast RT-qPCR Master (2x conc.)

• **Melting Analysis (FMCA)**

AddMelt qPCR Master with UDG (2x conc.)  
AddMelt RT-qPCR Master with UDG (2x conc.)

**Sample Preparation (spin column)** 15

• **Plasmid Extraction and Fragment Purification**

AddPrep Plasmid Extraction Kit  
AddPrep PCR Purification Kit  
AddPrep Gel Purification Kit  
AddPrep Fragment Purification Kit

• **Genomic DNA Extraction**

AddPrep Blood Genomic DNA Extraction Kit  
AddPrep Swab Genomic DNA Extraction Kit  
AddPrep Tissue Genomic DNA Extraction Kit  
AddPrep Plant Genomic DNA Extraction Kit  
AddPrep Bacterial Genomic DNA Extraction Kit  
AddPrep Soil Genomic DNA Extraction Kit

• **RNA Extraction**

AddPrep Total RNA Extraction Kit  
AddPrep Viral Nucleic Acid Extraction Kit

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
### Performance Data :

Product Name	Product Code	Provided with
Add Taq DNA Polymerase	17500	Add Taq DNA Polymerase (5 U/µl), 100 µl 10x Reaction Buffer 1.0 mL 10mM dNTP Mix. (2.5 mM each) 1.0 mL
Add Taq DNA Master (2x Conc.)	35001 (tube type)	1.0 mL x 5 tubes
	35001SS (8-strip type)	10.0 µl / tube x 96 tests

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

### Performance Data :

Product Name	Product Code	Provided with
AddMulti Taq DNA Master without UDG (2x conc.)	37101 (tube type)	1.0 mL x 5 tubes
	37101SS (8-strip type)	10.0 µl / tube x 96 tests
AddMulti Taq DNA Master with UDG (2x conc.)	37101U (tube type)	1.0 mL x 5 tubes
	37101USS (8-strip type)	10.0 µl / tube x 96 tests

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
### Performance Data :

Product Name	Product Code	Provided with
AddStart HP Taq DNA Polymerase	19101	AddStart HP Taq DNA Polymerase (5 U/µl), 100 µl 10x Reaction Buffer 1.0 mL 10mM dNTP Mix. (2.5 mM each) 1.0 mL
AddStart HP Taq DNA Master (2x conc.)	38101 (tube type)	1.0 mL x 5 tubes
	38101SS (8-strip type)	10.0 µl / tube x 96 tests

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## Reverse Transcription

- **cDNA Synthesis**  
AddScript cDNA Synthesis Kit
- **One-step RT-PCR**  
AddScript RT-PCR Master (2x conc.)

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
- Applications :**
- First - strand cDNA synthesis
  - cDNA library construction
  - Gene expression analysis

### Performance Data :

Product Name	Product Code	Provided with
AddScript cDNA Synthesis Kit	22701	20x Enzyme Solution 50 µl 5x Reaction Buffer 0.2 mL 10mM dNTP Mix. (2.5 mM each) 0.1 mL 10x Oligo dT20 0.1 mL 10x Random Hexamer 0.1 mL Nuclease-free Water 1.0 mL

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

## AddScript RT-PCR Master (without or with UDG), 2x conc.

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### Performance Data :

Product Name	Product Code	Provided with
AddScript RT-PCR Master without UDG (2x Conc.)	24701 (tube type)	1.0 mL x 5 tubes
	24701S (8-strip type)	10.0 µl / tube x 96 tests
AddScript RT-PCR Master with UDG (2x conc.)	24701U (tube type)	1.0 mL x 5 tubes
	24701US (8-strip type)	10.0 µl / tube x 96 tests

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

#### Performance Data :

Real-time PCR with bacterial gDNA (10<sup>5</sup>, 10<sup>4</sup>, 10<sup>3</sup>, 10<sup>2</sup> and 10 copies) by AddGreen qPCR Master

Real-time one-step RT-qPCR with human total RNA(2ng, 200pg, 20pg and 2pg) and human housekeeping gene (GAPDH(-) and ACTB(O) by AddGreen RT-qPCR Master

Product Name	Product Code	Provided with
AddGreen qPCR Master (2x conc.)	72501	1.0 mL x 1 tube
AddGreen RT-qPCR Master (2x conc.)	71302	1.0 mL x 1 tube

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

#### Performance Data :

Multi-plex Real-time PCR with bacterial gDNA (5x10<sup>5</sup>, 5x10<sup>4</sup>, 5x10<sup>3</sup>, 5x10<sup>2</sup> and 50 copies) by AddProbe qPCR Master with UDG (Cat. No. 72401U)

Quantitative multi-plex Real-time RT-qPCR with viral RNA (10<sup>6</sup>, 10<sup>5</sup>, 10<sup>4</sup>, 10<sup>3</sup>, 10<sup>2</sup> and 10 copies) by AddProbe RT-qPCR Master with UDG (Cat. No. 74201U)

Product Name	Product Code	Provided with
AddProbe qPCR Master (2x conc.)	72401	1.0 mL x 1 tube
AddProbe qPCR Master with UDG (2x conc.)	72401U	1.0 mL x 1 tube
AddProbe RT-qPCR Master (2x conc.)	74201	1.0 mL x 1 tube
AddProbe RT-qPCR Master with UDG (2x conc.)	74201U	1.0 mL x 1 tube

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### AddFast qPCR Master (without or with UDG), 2x conc.

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
### AddFast RT-qPCR Master (without or with UDG), 2x conc.

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
#### Performance Data :

Comparison test of AddProbe qPCR Master with UDG (Cat. No. 72401U) and AddFast qPCR Master with UDG (Cat. No.72401TNU) with human gDNA (20ng, 2ng, 200pg, 20pg) by multi-plex Real-time PCR

- PCR condition of AddProbe qPCR Master: 95°C 10 min, 95 °C 10 sec → 60 °C 30sec (40 cycles), total running time: 1hr 10min
- PCR condition of AddFast qPCR Master: 95°C 3 min, 95 °C 1 sec → 60 °C 3sec (40 cycles), total running time: 40min

Comparison test of AddProbe RT-qPCR Master with UDG (Cat. No. 74201U) and AddFast RT-qPCR Master with UDG (Cat. No.74201TNU) with viral RNA 5x10<sup>5</sup>, 5x10<sup>4</sup>, 5x10<sup>3</sup>, 5x10<sup>2</sup> and 50 copies) by multi-plex Real-time PCR

- Fluorescence dye: blue line – FAM, green line – HEX, red line: Texas Red
- RT-PCR condition of AddProbe RT-qPCR Master: 50 °C 20 min, 95°C 10 min, 95 °C 25 sec → 60 °C 50sec (40 cycles), total running time: 1hr 48min
- RT-PCR condition of AddFast RT-qPCR Master: 55 °C 5 min, 95°C 3 min, 95 °C 5 sec → 60 °C 10sec (40 cycles), total running time: 1hr 8min

Product Name	Product Code	Provided with
AddFast qPCR Master (2x conc.)	72401TN	1.0 mL x 1 tube
AddFast qPCR Master with UDG (2x conc.)	72401TNU	1.0 mL x 1 tube
AddFast RT-qPCR Master (2x conc.)	74201TN	1.0 mL x 1 tube
AddFast RT-qPCR Master with UDG (2x conc.)	74201TNU	1.0 mL x 1 tube

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
### Performance Data :

Comparison test of Fluorescence Melting Curve Analysis (FMCA) with human gDNA and two primer/probe(HEX) set for 2 targets detection between Normal qPCR Master (Competitor, blue line) and AddMelt qPCR Master(Addbiomeditek Co., red circle line)

Product Name	Product Code	Provided with
AddMelt qPCR Master with UDG (2x conc.)	72401M	1.0 mL x 1 tube
AddMelt RT-qPCR Master with UDG (2x conc.)	74201M	1.0 mL x 1 tube

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## Sample Preparation (Spin column)

- **Plasmid Extraction and Fragment Purification**

- AddPrep Plasmid Extraction Kit
- AddPrep PCR Purification Kit
- AddPrep Gel Purification Kit
- AddPrep Fragment Purification Kit

- **Genomic DNA Extraction**

- AddPrep Blood Genomic DNA Extraction Kit
- AddPrep Swab Genomic DNA Extraction Kit
- AddPrep Tissue Genomic DNA Extraction Kit
- AddPrep Plant Genomic DNA Extraction Kit
- AddPrep Bacterial Genomic DNA Extraction Kit
- AddPrep Soil Genomic DNA Extraction Kit

- **RNA Extraction**

- AddPrep Total RNA Extraction Kit
- AddPrep Viral Nucleic Acid Extraction Kit

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#### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Washing	15 ml
Binding	50 ml	Elution	20 ml

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

### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Washing	15 ml
Binding	60 ml	Elution	20 ml


### Performance Data :

Plasmids was extracted from 3mL overnight cultured recombinant E. coli by AddPrep Plasmid DNA Extraction Kit and loading volume is 2 µl extracted plasmid from 100 µl eluates.

The 20 µl PCR products which various amplicon size were purified by AddPrep PCR Purification Kit and loading volume is 5 µl purified PCR product from 50 µl eluates.

The 20 µl PCR products which various amplicon size were purified by AddPrep Fragment Purification Kit and loading volume is 5 µl purified PCR product from 50 µl eluates.

Product Name	Product Code	Provided with
AddPrep Plasmid Extraction Kit	10010	100 prep.
AddPrep PCR Purification Kit	10078	100 prep.
AddPrep Gel Purification Kit	10098	100 prep.
AddPrep Fragment Purification Kit	10088	100 prep.

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### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Binding	25 ml
Washing 1	30 ml	Washing 2	12 ml
Elution	25 ml	Proteinase K (20 mg/ml)	1.2 ml X 2 tubes

## AddPrep Swab Genomic DNA Extraction Kit, 100 prep.

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### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Lysis	30 ml
Binding	25 ml	Washing 1	30 ml
Washing 2	12 ml	Elution	25 ml
Proteinase K (20 mg/ml)	1.2 ml X 2 tubes		

## AddPrep Tissue Genomic DNA Extraction Kit, 100 prep.

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### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Lysis	30 ml
Binding	25 ml	Washing 1	30 ml
Washing 2	12 ml	Elution	25 ml
Proteinase K (20 mg/ml)	1.2 ml X 2 tubes		

## AddPrep Plant Genomic DNA Extraction Kit, 100 prep.

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### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Lysis	30 ml
Precipitation	20 ml	Binding	25 ml
Washing 1	30 ml	Washing 2	12 ml
Elution	25 ml	Proteinase K (20 mg/ml)	1.2 ml X 2 tubes

## AddPrep Bacterial Genomic DNA Extraction Kit, 100 prep.

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### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Lysis	25 ml
Binding	25 ml	Washing 1	30 ml
Washing 2	12 ml	Elution	25 ml
Lysozyme Buffer	55 ml	Lysozyme Solution (50 mg/ml)	1.2 ml X 2 tubes
Proteinase K (20 mg/ml)	1.2 ml X 2 tubes		

## AddPrep Soil Genomic DNA Extraction Kit, 100 prep.

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


### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Lysis	55 ml
Precipitation	22 ml	Binding	32 ml
Washing 1	30 ml	Washing 2	12 ml
Elution	25 ml	Lysozyme Buffer	55 ml
Lysozyme Solution (50 mg/ml)	1.2 ml X 2 tubes	Proteinase K (20 mg/ml)	1.2 ml X 2 tubes

**Performance Data :**

Total DNA was purified from 100mg soil samples by three different soil DNA Extraction Kits and 5 µl extracted gDNA from 100 µl eluates.

- Lane A: Competitor A
- Lane B: Competitor B
- Lane C: AddPrep Soil Genomic DNA Extraction Kit

Product Name	Product Code	Provided with
AddPrep Blood Genomic DNA Extraction Kit	10023B	100 prep.
AddPrep Swab Genomic DNA Extraction Kit	10023BS	100 prep.
AddPrep Tissue Genomic DNA Extraction Kit	10023T	100 prep.
AddPrep Plant Genomic DNA Extraction Kit	10023P	100 prep.
AddPrep Bacterial Genomic DNA Extraction Kit	10027	100 prep.
AddPrep Soil Genomic DNA Extraction Kit	10056	100 prep.

## AddPrep Total RNA Extraction Kit, 50 prep.

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### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column 1 (White ring)	50 ea	Spin column 2 (Green ring)	50 ea
Blood Lysis	75 ml	Lysis	25 ml
Binding	50 ml	Washing 1	30 ml
Washing 2	10 ml	Elution	25 ml
DNase I Reaction Buffer	1.0 ml x 3 tubes	Proteinase K (20 mg/ml)	1.2 ml x 1 tube
DNase I (1 U/µl)	0.5 ml x 1 tube		

## AddPrep Viral Nucleic Acid Extraction Kit, 100 prep.

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
### Kit Components

Solution & Material	Size	Solution & Material	Size
Spin column	100 ea	Lysis	40 ml
Washing 1	30 ml	Washing 2	12 ml
Elution	20 ml		

Human total RNA was purified from 200 µl whole blood sample by two different Total RNA Extraction Kits. A one-step RT-PCR was performed with 1 µl of 10-fold diluted eluate as RT-PCR template, various primers using AddScript RT-PCR Master (Addbiomeditek. Co.)

Lane A: human ACTB (505bp, 890bp)

Lane B: human HPRT (502bp, 863bp)

Lane C: human TUBB (450bp, 743bp)

Lane D: human IL (460bp, 535bp)

Product Name	Product Code	Provided with
AddPrep Total RNA Extraction Kit	10119	50 prep.
AddPrep Viral Nucleic Acid Extraction Kit	10034	100 prep.

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BEYOND BIOTECHNOLOGY  
BEYOND INNOVATION

CONTACT US

FOR BETTER HEALTH CARE  
**addbio**  
MEDITEK (주)애드바이오메디텍

대전광역시 유성구 테크노1로 41  
T. 042)934-2026 / F. 042)934-2029  
E. [addbiomeditek@addbiomeditek.com](mailto:addbiomeditek@addbiomeditek.com)